


instruction manual

Rodent Treadmill NG Cat. No. 47300



UGO BASILE S.R.L. sales@ugobasile.com / service@ugobasile.com <u>www.ugobasile.com</u>

UGO BASILE S.R.L.

Via Di Vittorio, 2 21036 GEMONIO, VA, ITALY Phone: +39 0332 744574

sales@ugobasile.com / service@ugobasile.com www.ugobasile.com

instruction manual

Rodent Treadmill NG Cat. No. 47300

Serial No.

SAFETY CONSIDERATIONS

ALTHOUGH THIS INSTRUMENT HAS BEEN DESIGNED WITH INTERNATIONAL SAFE-TY STANDARD, THIS MANUAL CONTAINS INFORMATION, CAUTIONS AND WARN-INGS WHICH MUST BE FOLLOWED TO ENSURE SAFE OPERATION AND TO RETAIN THE INSTRUMENT IN SAFE CONDITIONS.

SERVICE AND ADJUSTMENTS SHOULD BE CARRIED OUT BY QUALIFIED PERSONNEL, AUTHORIZED BY UGO BASILE ORGANIZATION.

ANY ADJUSTMENT, MAINTENANCE AND REPAIR OF THE OPENED INSTRUMENT UN-DER VOLTAGE SHOULD BE AVOIDED AS MUCH AS POSSIBLE AND, WHEN INEVITA-BLE, SHOULD BE CARRIED OUT BY A SKILLED PERSON WHO IS AWARE OF THE HAZARD INVOLVED.

CAPACITORS INSIDE THE INSTRUMENT MAY STILL BE CHARGED EVEN IF THE IN-STRUMENT HAS BEEN DISCONNECTED FROM ITS SOURCE OF SUPPLY.



www.ugobasile.com

Rodent Treadmill

Cat. No. 47302 for Rats Cat. No. 47303 for Mice



General

"Exercise is a multifactorial activity that affects virtually every organ and tissue in the body. Not only does exercise contribute many health benefits, but lack of exercise is implicated in many chronic health problems.


As evidence continues to accumulate concerning the impressive range of health benefits that exercise confers, biomedical researchers have increasingly become interested in conducting systematic studies of exercise to further define those benefits"

(from Resource Book for the Design of Animal Exercise Protocols, APS, Feb 2006)


Ugo Basile introduces an original TREADMILL for rats and mice. The same device is suitable for tests on either rats or mice, by simply replacing the lane assembly.

Our model incorporates a shock grid at the back of the treadmill to deliver a mild electric shock, when an aversive stimulus is required.

The running-lane assembly can be manually tilted from -25° to $+25^{\circ}$, in steps of 5° .



- MEASURES ENDURANCE, DISTANCE, SPEED
- SAME DEVICE TO TEST RATS & MICE
- COMPACT AND USER-FRIENDLY: test settings & monitoring controlled by the attached electronics and managed on the touch-screen.



Main Features

- SPEED: from 3 to 100m/min, in steps of 1m/min
- MODES: constant, accelerating, custom ramps
- SLOPE: positive (uphill) or negative (downhill), from -25° to +25°
- SHOCK: from 0 to 2mA (in 0.1mA steps), included
- CONTROLS: 4"3 touch-screen to set and monitor the test
- **X-PAD** SOFTWARE: brand new, user-friendly version, to set the experiment and manage the results
- DETECTION: via incorporated electronic circuit automatically detects speed & absolute and relative distances

Ugo Basile: more than 10,000 citations



CHECK-LIST

□ 47302 Rat Treadmill

Mouse Treadmill □ 47303

□ 47300 Combo Package Rat&Mouse

CLIENTE / CUSTOMER_____

Ordine No. / Order No. _____ Data / Date___ /____/

UB code	CAT.No.	47302	47303	47300	DESCRIPTION		DESCRIZIONE		V
	47301	1	1	1	MAIN UNIT: TAPIS	S –ROULANT	UNITÀ BASE TA	PIS-ROULANT	
	47300-002	1		1	RAT LANE ASSEMBLY WITH		ASSIEME CORS		
	47300-003		1	1			ASSIEME CORS		
E-WP 008		1	1	1	MAINS CABLE	EUROPE	CAVO RETE	EUROPA	
E-WP 008-1		I	I	I	MAINS CADLE	U.S.A.	CAVUREIE	U.S.A.	
E-FT 021		1	1	1	SET OF 2 FUSES	(T5A 5X20)	SET DI 2 FUSIBI	LI (T5A 5X20)	
E-AU 041 USB	47300-302				MANUAL		MANUALE		
pen-drive	X-PAD				SOFTWARE		SOFTWARE		
OPTIONAL	OPTIONAL								
	47300-013				6-LANE ASSEMBLY, WITHOUT LID, FOR TETHERED MICE)		ASSIEME 6 COF CABLATI, SENZ		
	1				1		1		

DATE	DATE / /		Serial No.	IMBALLATO DA / PACKED BY						
Set for	Set for									
	115V 5	0 Hz	230V 50 Hz	□ 115V 60 Hz □ 230V 60 Hz						

IMPORTANT/IMPORTANTE:							
Check the shipment for completeness immediately after receipt: should you find any discrepancy, please fill in the following part and transmit it to our fax no. +39 0332 745488							
Al ricevimento della merce controllate che la spedizione sia completa: in caso di discrepanza, completate il formulario di seguito riportato ed in- viatelo al nostro fax no. 0332 745488							
FROM: Name	Company/Institution						
DATE	REF.						
NOTE							
MOD.04 REV 0							



CONTENTS

1	GENERAL1	
1.1	PRINCIPLE OF OPERATION1	
2	INSTRUMENT DESCRIPTION	:
2.1 2.2 2.3 2.4 2.5 2.5.1 2.6 2.7 2.8	RUNNING BELT2RAT LANE-ASSEMBLY2MOUSE LANE-ASSEMBLY3SHOCKER CIRCUIT AND SHOCK GRID3CONTROLLER WITH TOUCH-SCREEN3NOTES ON RESISTIVE TOUCH-SCREENS4EXPERIMENTAL CONFIGURATION5DATA COLLECTION AND MANAGEMENT5TREADMILL FEATURES5	
3	INSTALLATION6	;
3.1 3.2 3.3 3.4 3.5 3.5.1 3.5.2 3.5.3 3.6 3.7 3.8 3.8.1 3.8.2 3.9	UNPACKING & PRELIMINARY CHECK6NOTES ON THE INSTRUCTION MANUAL6GENERAL SAFETY INSTRUCTIONS7ASSEMBLING THE INSTRUMENT7BEFORE APPLYING POWER7MAINS SWITCH7FUSE HOLDER7MAINS CORD7INTENDED USE7ADDITIONAL SAFETY CONSIDERATION7CONNECTIONS9USB HARD DISK9INSTALLATION OF X-PAD SOFTWARE9	\$, , , , , , , , , , , , ,
4	PRELIMINARY9	1
4.1 4.2	PLACING THE ANIMAL ON THE RUNNING SURFACE	
5	OPERATION)
5.1 5.2 5.2.1 5.2.2 5.3 5.4	ANIMAL SELECTION. 10 HOME PAGE 11 SYSTEM ICONS 11 HOME PAGE ICONS 11 EXPERIMENT. 12 SET-UP 13	
5.4.1	Settings 13	ŝ



5.4.2	SPEED	13
5.4.3	CONSTANT SPEED	13
5.4.4	ACCELERATED/DECELERATED SPEED	14
5.4.5	CUSTOM RAMP	14
5.4.6	Shock	14
5.5	START BUTTON	15
5.5.1	STARTING AT CONSTANT SPEED	15
5.5.2	STARTING A RAMP	16
5.6	RESULTS	17
5.6.1	How to Export, Interpret and Erase the Experimental Data	17
5.7	IMPORTANT NOTE ON SHOCK FREQUENCY/DETECTION	18
5.8	USB STORAGE	19
5.8.1	SAVE DATA	19
5.8.2	Load Protocol	19
5.9	UTILITIES	20
5.9.1	Software Update	20
5.9.2	SETTING DATE & TIME	21
5.9.3	DIAGNOSTICS	22
5.9.4	Reset HD	22
5.9.5	CLEANING	23
5.10	SIGNAL OUTPUT	23
5.10.1	PIN-OUT	23
•		~ 4
6	MAINTENANCE	24
_		
6.1	ELECTRICAL	24
6.1 6.2	ELECTRICAL	24 24
6.1 6.2 6.3	ELECTRICAL CLEANING RUNNING BELT	24 24 24
6.1 6.2 6.3 6.3.1	ELECTRICAL. CLEANING RUNNING BELT. ALIGNING THE BELT	24 24 24 24
6.1 6.2 6.3 6.3.1 6.3.2	ELECTRICAL CLEANING RUNNING BELT ALIGNING THE BELT REGULATING THE BELT'S STRAIN	24 24 24 24 25
6.1 6.2 6.3 6.3.1 6.3.2 6.3.3	ELECTRICAL. CLEANING. RUNNING BELT. ALIGNING THE BELT REGULATING THE BELT'S STRAIN LANE ASSEMBLY	24 24 24 24 25 26
6.1 6.2 6.3 6.3.1 6.3.2 6.3.3 6.3.4	ELECTRICAL. CLEANING RUNNING BELT. ALIGNING THE BELT REGULATING THE BELT'S STRAIN LANE ASSEMBLY DROP PAN	24 24 24 25 26 27
6.1 6.2 6.3 6.3.1 6.3.2 6.3.3 6.3.4 6.3.5	ELECTRICAL. CLEANING RUNNING BELT. ALIGNING THE BELT REGULATING THE BELT'S STRAIN LANE ASSEMBLY DROP PAN GRID ASSEMBLY	24 24 24 25 26 27 28
6.1 6.2 6.3 6.3.1 6.3.2 6.3.3 6.3.4 6.3.5 6.4	ELECTRICAL CLEANING RUNNING BELT ALIGNING THE BELT REGULATING THE BELT'S STRAIN LANE ASSEMBLY DROP PAN GRID ASSEMBLY LUBRICATION	24 24 24 25 26 27 28 29
6.1 6.2 6.3 6.3.1 6.3.2 6.3.3 6.3.4 6.3.5 6.4 6.5	ELECTRICAL. CLEANING RUNNING BELT ALIGNING THE BELT ALIGNING THE BELT REGULATING THE BELT'S STRAIN LANE ASSEMBLY DROP PAN GRID ASSEMBLY LUBRICATION LONG INACTIVITY	24 24 24 25 26 27 28 29 29
6.1 6.2 6.3 6.3.1 6.3.2 6.3.3 6.3.4 6.3.5 6.4	ELECTRICAL CLEANING RUNNING BELT ALIGNING THE BELT REGULATING THE BELT'S STRAIN LANE ASSEMBLY DROP PAN GRID ASSEMBLY LUBRICATION	24 24 24 25 26 27 28 29 29
6.1 6.2 6.3 6.3.1 6.3.2 6.3.3 6.3.4 6.3.5 6.4 6.5 6.6	ELECTRICAL CLEANING RUNNING BELT ALIGNING THE BELT ALIGNING THE BELT REGULATING THE BELT'S STRAIN LANE ASSEMBLY DROP PAN GRID ASSEMBLY LUBRICATION LONG INACTIVITY CUSTOMER SUPPORT	24 24 25 26 27 28 29 29 29
6.1 6.2 6.3 6.3.1 6.3.2 6.3.3 6.3.4 6.3.5 6.4 6.5	ELECTRICAL. CLEANING RUNNING BELT ALIGNING THE BELT ALIGNING THE BELT REGULATING THE BELT'S STRAIN LANE ASSEMBLY DROP PAN GRID ASSEMBLY LUBRICATION LONG INACTIVITY	24 24 25 26 27 28 29 29 29
6.1 6.2 6.3 6.3.1 6.3.2 6.3.3 6.3.4 6.3.5 6.4 6.5 6.6 7	ELECTRICAL CLEANING RUNNING BELT ALIGNING THE BELT REGULATING THE BELT'S STRAIN LANE ASSEMBLY DROP PAN GRID ASSEMBLY LUBRICATION LONG INACTIVITY CUSTOMER SUPPORT 47300 SPECIFICATIONS	24 24 24 25 26 27 28 29 29 29 29 30
6.1 6.2 6.3 6.3.1 6.3.2 6.3.3 6.3.4 6.3.5 6.4 6.5 6.6	ELECTRICAL CLEANING RUNNING BELT ALIGNING THE BELT ALIGNING THE BELT REGULATING THE BELT'S STRAIN LANE ASSEMBLY DROP PAN GRID ASSEMBLY LUBRICATION LONG INACTIVITY CUSTOMER SUPPORT	24 24 24 25 26 27 28 29 29 29 29 30
6.1 6.2 6.3 6.3.1 6.3.2 6.3.3 6.3.4 6.3.5 6.4 6.5 6.6 7	ELECTRICAL CLEANING RUNNING BELT ALIGNING THE BELT REGULATING THE BELT'S STRAIN LANE ASSEMBLY DROP PAN GRID ASSEMBLY LUBRICATION LONG INACTIVITY CUSTOMER SUPPORT 47300 SPECIFICATIONS	24 24 24 25 26 27 28 29 29 29 29 30 31
6.1 6.2 6.3 6.3.1 6.3.2 6.3.3 6.3.4 6.3.5 6.4 6.5 6.6 7	ELECTRICAL. CLEANING RUNNING BELT ALIGNING THE BELT ALIGNING THE BELT REGULATING THE BELT'S STRAIN LANE ASSEMBLY DROP PAN GRID ASSEMBLY LUBRICATION LONG INACTIVITY CUSTOMER SUPPORT 47300 SPECIFICATIONS ORDERING INFORMATION	24 24 24 25 26 27 28 29 29 29 29 30 31

FIGURE INDEX



Rodent Treadmill NG

Cat. 47300

1 GENERAL

"Exercise is a multifactorial activity that affects virtually every organ and tissue in the body. Not only does exercise contribute many health benefits, but lack of exercise is implicated in many chronic health problems.

As evidence continues to accumulate concerning the impressive range of health benefits that exercise confers, biomedical researchers have increasingly become interested in conducting systematic studies of exercise to further define those benefits".¹

Fatigue is a common and frequently poorly-understood symptom in many diseases and disorders. New preclinical assays of fatigue may help to improve current understanding of fatigue-like behavior in rodents and many other exercise paradigms and study future treatment of fatigue.

Ugo Basile introduces an original TREADMILL for rodents: the same device is suitable for tests on either rats or mice, by simply replacing the lane assembly.

The 47300 Treadmill is a compact and user-friendly device: test settings & monitoring are managed by the attached control unit. It measures endurance, distance (absolute and relative) and speed.

1.1 Principle of Operation

Treadmills are rolling belts (tapis-roulants) with presettable speed and adjustable uphill and downhill inclination (slope), enabling forced exercise training and accurate testing of fatigue in lab animals.

"Treadmill running has been used extensively over the past decades to study behavioral, physiological, biochemical, and, more recently, molecular responses to both acute exercise stress and chronic exercise training. Although investigators have used a wide variety of species (...) for treadmill running studies, they have used rodents in most of these studies.

Treadmill running has the distinct advantage over other forms of exercise, including spontaneous wheel running and swimming, that the total amount of external work done by the rat can be easily calculated⁷¹

¹ from Resource Book for the Design of Animal Exercise Protocols, APS, Feb 2006



2 INSTRUMENT DESCRIPTION

Our Treadmill consists of a main unit, incorporating the drive, the shocker, the running belt and the control unit with 4"3 touch-screen.

The assembly also incorporates a shock grid, positioned at the bottom of the ramp, made of 3mm bars, spaced 8mm apart.

Two different lane assemblies are available, to provide the ideal running tracks for either rats or mice: replacing the lane assembly is extremely easy.

The whole hardware, except the lane assembly, is the same for both rats and mice: animal selection is quickly done when switching on the device, see paragraph 5.1-Animal Selection.

The running surface can be manually tilted from -25° to +25°.



Figure 1 "Mouse Lane-Assembly"

The device features an auto-cleaning tool, see paragraph 5.9.5, and a pan to collect droppings, see paragraph 6.3.4.

Our treadmills are designed for optimal performance under conditions of intensive use.

2.1 Running Belt

Rodents running on the treadmill must be able to maintain good traction while walking or running, to prevent slipping. The running surface of the 47300 consists of an easy-to-clean alimentary-grade white belt, providing suitable grip.

The walking surface is not porous and is soft enough to eliminate toenail and foot problems that may arise from daily bouts of exercise training.

This specially selected material make the treadmills easy to clean and disinfect and require minimal maintenance, see paragraphs 6.2-Cleaning.

2.2 Rat Lane-Assembly

The rat assembly, a structure which is quickly and easily fitted to the main unit, consists of 15cm high external walls and inside partitions, to divide the running belt into 3 lanes, each 45x11cm. Each lane is provided with a transparent lid.

Removing and replacing the lane assembly is an easy job: see paragraph 6.3.1.



2.3 Mouse Lane-Assembly

The Mouse assembly, similar to the rat device, has different dimensions: walls and partitions are 7cm high, and the running belt is divided into 6 lanes, each 45x5.5cm. Each lane is provided with a transparent lid.

A special lane assembly for tethered mice (47300-013) is also available: each lane 45x5.5, height 15cm, without lid, see paragraph 8.1.

2.4 Shocker Circuit and Shock Grid

"Treadmill running may be construed as a form of forced exercise in which the animal does not have a choice of participating in the activity. Because of this, noxious stimuli (e.g., electric shock and bursts of high-pressure air) may be needed to motivate the animals to exercise" ¹

Our model incorporates as standard a shocker to deliver a mild electric shock, when an aversive stimulus is required, via deliver the shock grid positioned at the bottom of the walking slope.

The shock circuit is incorporated in the main unit. Shock intensity and frequency can be preset via the attached controller module (from 0 to 2mA, 1, 2 or 3Hz), as well as the cut-off number of shocks. The setting is common to all lanes.

The grid attached at the bottom of the running belt delivers the foot-shock and also functions as detection circuit.

The assembly incorporating the grid, its circuit, and the array of START buttons can be easily removed for cleaning: see paragraph 6.3.4.

The same grid also functions as detection system: the controller detects absolute and relative distances, which are recorded together with speed.



Figure 2 "Grid Assembly"


2.5 Controller with Touch-Screen

The controller is part of the main unit, attached to the tapis-roulant.

It incorporates a 4"3 touch-screen, for basic setting and monitoring, via an intuitive panel.



The picture shows the main menu.



During the test, the touch-screen indicates the number of shocks received, the elapsed time, and the walked distance.

When the test is started the lane background colour, which is normally green, becomes light blue.

The test results for each section, can be browsed in a condensed view.



The brand new **X-PAD** software, included as standard, provides a user-friendly interface to set the experiment and a powerful tool to manage the results easily on one's PC. Data exchange from the instrument to the PC is done via the USB flash drive provided.

See also paragraphs 2.6-Experimental Configuration and 2.7-Data Collection and Management.

2.5.1 Notes on Resistive Touch-Screens

If you have in mind the sensitivity to touch of your tablet or smart-phone, you might think the controller touch screen does not react promptly, and you have to press down harder.



This is because we employ a resistive touch-screen (while most smart phones use a capacitive screen); the resistive screens are better of our application because of

- High resistance to dust and water
- Better use with gloved hand or stylus (because the conductive technology is dependent on the conductive nature of human body, it doesn't work if the user is wearing gloves).

2.6 Experimental Configuration

Via the new *X-PAD* software, the operator can easily organize the experiment on her/his PC, and upload it to the Treadmill via the USB flash drive provided.

Treatments, protocols, stages, animals, and various test features (speed, mode, revolution, etc.) can be quickly defined and saved for future use.

di d		See.					TOOL	BAR				
ana				1 Tests	-	-			200		_	_
Contract		Annet		Ever D	Minul ()	Anna Cone		inge ig beareast	Permit	Na Sing	Data - Terra	
Suprat 2	Tearrent Tea 3	(20				Josen .	Main 🗖			1.00	4416-001	0
Ceatries 1	Testret Tet 1	- 22				\$cont.	Mate 🖸	EDC Textnert1	0 0	3.00	Ga 04/05/2016	8
Teatrent 2	Testrert Test 7					(room)	Main 🖬			100	04/06/2019	
Contract I	TestretTet1	10				(beaution)	Maria 🖸		0 0	139	Cente-2016	
autorit 2						Brown .	Maie 🖸			100	Devision	0
Instruct S						Becom.	Maria 🗖		0 0	1.00	Dete 275	
internal internal	Teation fant			11 - 11		Boun	Maie 🖬		0 0	100	D4/10/2015	
ament					1 1	Desire	Xilia 🗖			1.00	Cu-06/27%	
enis Testa				1.1	4 1	parts.	Note 🖸		0 0	1.02	C 04/06/2275	0
Total State				2.3		Doctor 1	State 🖬	5.00 Testment 8	0 0	1.01	D 0498.075	-
					1.10	prises.	Main 🖬	0.00 Teatman 10	0 0	1.00	04/06/0015	
Real and the second	Tare					(trans	Mare D	520 Testnert 11	0 0	100	0415-275	
share 2						\$tpan	Main D	0.00 Testnert V	0 0	100	C 0408-2015	
						Sec. 1	Max 🖬		ā ā	120	C4/06/2275	
			Duality			Score .	the D		0 0	150	Dertrictors.	8
tope 2	2008 Per 2		10.00			2 mart	Main C		ă ă	1.00	0406005	ö
inge 2	Date Terr 1		10.00			Press.	Main C		0 0	1.00	CA-06.00%	8
inge Z	Date Test 2		10.00			and a second			8	10	CA16.001	
Sec.)	Dage Sect 2		2000			low l	Maia 🖬			100	0+16/225	No.
and a			-			Brown .	New O		ă ă	100	Detects	ă
				_					0 0	1.00		
	DOCK					Print .	Note 🖬				04/05/2019	
	JULK	AR		1.000		lation .	the			1.00	0448.20%	
			and the second second	1.000		bovi	Maine .				C+16.001	
				·		Scon	tine 🖬	4.00 74			C+06305	
			And in case of the local division in which the local division in which the local division in which the local division is not the local division in the loc	· · ·		Stown	Main 🖬	0.00 Te	ΛΛΙΝΙ	ADEA	Cetterrn	
47.		officers. Total T		1.1		Print.	Mate 🖬	10.74	//AUNI	AREA	Cuss.275	
nje Z	Dage Test D.D.		6.00	1.0		\$ dam	Mala 🖬				C4/06/2011	
ispe 1	Dege Text D-L		0.00	·	2 2	Brian .	Mate 🖬	01.74			Currences.	-
isge 2	Dage Text 2-2		0.00			\$pen	Mate 🖬	0.00 Technologian		- 04	G4/06/2015	
ige 2			000	1.00	27 3	Sec.	Male 🖬	200 Teament 29	0 0	100	D455075	
					20 3	Irpan	Male 🖸	000 Teatment 30	0 0	139	04/06/2015	
1111				1	M 34	press.	Main 🖸	U.O. Textment H	0 0	1.00	Col 04 05 27%	
			Deather		34 11	(inclusion)	Male D	0.00 Teatment 32	0 0	100	04/06/2016	
ige 0			5.00			print.	Male 🖸	100 Testnert II	0 0	100	Di-06/2015	10
tige 1					14 14	Series .	Male 🖸	0.00 Testment M	0 0	100	C4/06/2015	8
and a second sec					- 25 - 24	Sec	Male 🖸	0.00 Seatment SI.	0 0	1.50	C4/06/2016	8
							his D		0 0	1.00	04/06/2015	8
						Print I	1000			100	Contraction of the	ă i

w

w

X

Configurations are easily:

- exported to Text, Excel or Pdf reports
- saved to cloud to DropBox, OneDrive, GoogleDrive State

2.7 Data Collection and Management

A basic version of the collected data can be viewed on the touch-screen; when transferred to PC via USB drive, test results appear in full version.

The software automatically classifies the data, combining configuration settings with test results. The user can add information, before or after the test.

Results appear in a tree-like structure, where columns can be dragged and dropped to customize the layout.

Data are quickly:

- exported to Text, Excel or Pdf reports
- saved to cloud to DropBox, OneDrive, GoogleDrive State

2.8 Treadmill Features

The main features of the new model are:

• **SPEED**: adjustable in the range 3-100m/min, in steps of 1 m/min



- MODE: constant, ramp (accelerating), multi-step ramp
- **SHOCK**: 0 to 2mA, 1, 2 or 3Hz
- **SLOPE**: uphill or downhill, manual tilting, from -25° to +25°, in steps of 5°
- **CONTROLS**: 4"3 touch-screen to set and monitor the test
- SOFTWARE: X-PAD brand new, user-friendly version, to set the experiment and manage the results
- **DETECTION**: via grid at the bottom of the running belt

3 INSTALLATION

3.1 Unpacking & Preliminary Check

Check the contents of the shipment for completeness, packing list to hand, and visually inspect the instrument as soon you take it out of the packaging. Use the supplied *Check List*.

If the instrument is damaged, inform the carrier immediately, notifying our company. If after having tested it, the instrument fails to meet rated performances, please contact our after sales service, see paragraph 6.6-Customer Support.



The 47300 is shipped in a dedicated wooden crate: we suggest to store the crate and the dedicated packaging material for stocking and possible future shipment.



Protect the environment!

If necessary, dispose of packaging properly, according to existing and applicable waste management rules and regulations.

3.2 Notes on the Instruction Manual

The Instruction Manual included in the package (on the USB pen drive) is necessary for the correct installation and operation of the instrument.

We recommend reading the manual with attention, as it is essential for the correct installation and operation of the instrument.

Please save the manual, ready to be consulted by the qualified personnel who use the instrument. Print it, only if necessary.

Our Instruction Manuals are available as free download on our web. For any additional information and/or assistance, you are welcome to contact our Service Department (see paragraph 6.6-Customer Support), specifying the serial number of your instrument.



3.3 General Safety Instructions

The following guidelines must be followed to ensure safe operation.

- **! DO NOT** attempt to open or perform any service work
- ! DO NOT connect up human subjects

3.4 Assembling the Instrument

The Treadmill is shipped assembled: position it on a stable, sturdy, and reasonably flat bench or table surface. It is a heavy device: please remember this when lifting it!

3.5 Before Applying Power

Consider the Power Module (see figure 3), positioned left on the back panel, which encompasses – from left to right - the mains switch, the inlet connection of the mains cord, and the fuse holder/voltage selector.

Figure 4, "Power Module"



3.5.1 Mains Switch

This two-pole toggle switch, which complies with international safety standards, is lighted when the instrument is **ON**.

3.5.2 Fuse Holder

The fuse holder comprises two fuses, one on the live, and the other on the neutral. Use (T5A) timed fuses for operation at both 115 or 230 Volts.

For fuse replacement, please refer to paragraph 6.1-Electrical.

3.5.3 Mains Cord

It is a standard cable, Cat. # E-WP008. Make sure your power outtake is provided with a reliable ground connection.

3.6 Intended Use

The Treadmill is intended for investigation use on **laboratory animals only**.

3.7 Additional Safety Consideration

1) Use original accessories and spare parts only, see also paragraph 8-ORDERING INFORMATION.





- 2) immediately disconnect and replace damaged mains cord.
- 3) do not obstruct a comfortable access to the power module.
- 4) do not operate in hazardous environments or outside prescribed environmental limitations (i.e. +10c° / +40c°, 95% max. relative humidity, non-condensing), see also paragraph 7-47300 SPECIFICATIONS.
- 5) do not spray any liquid on the connectors and on the geared motor.

UGO BASILE DOES NOT ACCEPT ANY RESPONSIBILITY FOR PROBLEMS OR HARM CAUSED TO THINGS OR PERSONS, ARISING FROM:

- incorrect electrical supply;
- incorrect installation procedure;
- incorrect or improper use or, in any case, not in accordance with the purpose for which the instrument has been designed and the warnings stated in the instruction manual supplied with the instrument;
- replacement of original components, accessories or parts with others not approved by the manufacturer;
- servicing carried out by unauthorized personnel.

3.8 Connections

Connect the mains cord to a power outtake, provided with a reliable earth connection.

Now consider the connection module: for easy and quick access, it is positioned on the front panel.

It encompasses the following connectors, from left to right:




Figure 5 "Connection Module"

- the upper USB enables data exchange (configurations & results) with the PC, and allows firmware upgrades, via the USB flash drive provided, or to connect a keyboard (see paragraph 3.8.2)
- the lower USB accommodates the USB storage key. WARNING: DO NOT <u>REMOVE IT! See paragraph 3.8.1</u>
- I/O this D-SUB 15 connector provides TTL outputs for shock status for each lane, and speed, see paragraph 5.10.
- COM this USB-B 2.0 allows communication to the PC (for Ugo Basile use, programming and debug)
- **ETH** the Ethernet connector will be used for remote diagnosis and Internet access used for standard operation.



3.8.1 USB Hard Disk



The USB storage key must **ALWAYS** be connected to the bottom port (see picture below). This USB key works as the inner memory of the system. It contains the collected data and the protocols.

If the device does not detect the USB HARD DISK, a message box asks to disconnect and to reconnect the USB Key and restart the device.



If during a normal operation, the user disconnects the USB HARD DISK, the device automatically switches off and restarts.

WARNING: this operation could damage the experiment data and impair the internal memory.

3.8.2 Keyboard

When setting up the experiment, it is often requested (although not compulsory) to enter data via the touch-screen. For example information about the test, stage, animal, treatment, etc.

To ease the data input, it is possible to connect a <u>standard QWERTY keyboard</u> to the upper USB port. The keyboard will be automatically recognized and **alphanumeric** data typed on the keyboard will appear on the touch-screen.

Beside all the alphanumeric keys, the PC will accept:

- The **backspace**, to cancel the entered data
- The **ESC** key, to leave the menu
- The ENTER key, to confirm the input and move to the following field

3.9 Installation of X-PAD Software

The *X-PAD* software is saved on the USB Flash Drive provided with the instrument: download it on the Hard Disk of your PC.

4 **PRELIMINARY**

4.1 Placing the Animal on the Running Surface

We recommend that **the belt is set in motion before placing the mice in position**, or by the time the last mouse is in place, the first may well be facing the wrong direction.

In constant speed mode, start the motion at the selected speed and then place the mice/rats, one by one, in their respective lanes, at the same time starting the related counters to zero, by the corresponding button.



Experience and patience will teach the ideal technique: as in most behavioural devices, the man/animal interface is not less important than the animal/machine one, and requires some patience and gentle handling to be perfected.

Speaking about animal/machine interface, rodents running on the treadmill must be able to maintain good traction while walking or running, to prevent slipping. The running surface of the 47300 consists of an easy-to-clean alimentary-grade white belt, providing suitable grip; the walking surface is not porous and is soft enough to eliminate toenail and foot problems that may arise from daily bouts of exercise training.

This specially selected material make the treadmills easy to clean and disinfect, and require minimal maintenance, see also paragraph 6.2. This type of belts usually wear well, with minimal breakage, and rarely need to be replaced during the lifetime of the treadmill.

4.2 Acclimation

The Treadmill is basically an exercise machine, which does not require the animal to learn a specific task. A proper training phase is therefore not required.

However, it advisable to carry out an acclimation period with low and constant speed, in order for the naïve mouse to familiarize with the running belt: each mouse is placed into its section for a short period, before starting the actual test.

5 OPERATION


Switch on the Treadmill: the instrument performs a CPU check, which might take up to 60 seconds, during which the display will show:



5.1 Animal Selection

When the Treadmill is switched on, the first page that shows up is the one requiring the selection of the animals on test.

Since the whole hardware, except the lane assembly, is the same for both rats and mice, it's important to specify which animals the device is going to deal with.



After the selection, it will not be possible to go back to this page unless the Treadmill is switched off and then on again.



5.2 Home Page

The test is managed by the 4"3 touch-screen.

This is how the Home Page looks like.

From the Home Page, the logo button on the tool-bar, opens the "screen-saver" window





This window displays current date and time, and the version of the software on-board.

From any other menu, the logo button always brings back to the Home page.

5.2.1 System Icons



NEW SOFTWARE DETECTED: if this icon is present, a new version of SW has been detected on the USB Key.



ETHERNET DETECTED: an Ethernet connection has been detected

USB KEY DETECTED: if this icon is present, a USB has been detected on the USB port.

5.2.2 Home Page Icons



In the EXPERIMENT page, it is possible to enter information about the test, including treatment, protocol, stage, trial, see paragraph 5.3.



Speed and Mode are defined from the SET-UP page, see paragraph 5.4



goes to START Page, see paragraph 5.5





goes to RESULTS Page, see paragraph 5.6



goes to USB STORAGE Page, see paragraph 5.8



goes to UTILITIES Page, see paragraph 5.9

5.3 Experiment

Entering this section the user will find two different layouts depending on the animals that have been previously selected:

- 6 distinct tabs for the mice
- 3 tabs for the rats

In this phase the user can fill manually the experiment tags lane by lane (Treatment, Protocol, Stage, Inclination, Trial, ID).

Monday 9:41	y		-	(🛇 📑 🖥			
Field	Nom	e:	Field	l cont	ent				
Q	w	E	R	Т	- \	r l	J	ı c	P
12?!	Α	S	D	F	G	н	J	ĸ	L
BAC	к	z	x	С	v	в	N	М	ОК

	•	🚍 🔹	
TREATMENT			
PROTOCOL			
STAGE	I	NCL	
TRIAL		ID	
LANE 1 LAN	IE 2 LANE 3 LAN	IE 4 LANE 5 LA	NE 6

Monday 9:41	Y		ې 🖯						
Field	Name	: F	ield c	onter	nt				
1	2	3	4	5	6	7	8	9	0
AZ	@	:	;	-	-	#	1	1	&
BAC	к	•	, 1	•	!	6	*	%	ок

Enter a maximum of 12 characters in each field. TREATMENT, PROTOCOL, STAGE, INCLINATION and TRIAL are text boxes, where it's possible to enter an information.

All the information boxes are independent for each lane: select another lane and enter the related data. Repeat the same process for each lane.

In the RESULT table the entered information will be matched with the test results. Of course, the user may choose not to enter any information and leave the fields empty.



The Home Page icon brings back to the home page.





NOTE:

Whenever data are entered on the virtual keyboard, it is possible to use an external QWERTY keyboard instead, connected to the device's upper USB port, see also paragraph 3.8.2.

5.4 Set-Up

The parameters that characterize an experiment are set in this page: **Settings**, **Speed** (including **Mode**) and **Shock**.

Having selected the Speed and mode, press \bigcirc to go back to the Home page.

5.4.1 Settings

The user chooses the thresholds cut-off which lead to the end of a test, if the animal does not end it earlier.

In other words, the test is concluded when one of the following parameters is reached within a single trial:

- The number of shocks an animal can receive
- The maximum distance run
- The maximum time elapsed

Each value can be preset from 0 to 999.



5.4.2 Speed

In the second section, the selection among the speed mode and the appropriate set point(s) has to be done. The user may choose among the **CONSTANT**, **RAMP**, **CUS-TOM RAMP**, modes, by depressing the related button at the page bottom.

5.4.3 Constant Speed

With **Constant Speed** mode, the belt rotates CCW at a speed which remains constant, for the duration of the test.






Depress the icon to activate the numeric keyboard, and enter a speed between 3 and 100 m/min. Confirm by "OK".

5.4.4 Accelerated/Decelerated Speed

This mode allows testing the subjects with an **increasing/decreasing** speed. It is necessary to set the initial and final speed, and the ramp (the time the rotor takes to go from the initial to the final speed).





Depress the icon **I** to activate the numeric keyboard to enter:

- Initial speed: between 3 and 100 m/min
- final speed: between 3 and 100 m/min
 - if the initial speed is lower than the final speed, the speed will be accelerating: on the contrary a decelerated speed is obtained by setting an initial speed higher than the final one)
- ramp: from 1s to 999s (acceleration limits: maximum 100m/min)

5.4.5 Custom Ramp

This feature is not controlled on the touch-screen, and is to be set on the PC. The related button is not active. For instruction, please refer to the X-PAD manual.

5.4.6 Shock

The last section defines the shock stimulus which will be delivered during the trial, setting both the frequency (1/2/3Hz) and the intensity (0-2.0 mA).

The frequency will determine the number of shocks per second the animal will receive, 1, 2 or 3 shocks per second.





IMPORTANT:

The different shock frequency will also impact the capability to detect shortdistance shock events, see parag. 5.7-Important Note On Shock Frequency.




5.5 Start Button

5.5.1 Starting at Constant Speed

In the constant speed mode, during **IDLE** state, the motor and the counters are disabled.

Pressing the **START** button, activates the motor. The motor stops and the counters are disabled by pressing the **STOP** button.



As the walking surface is a single belt (not separate for each lane), when the motor is started all the lanes will be moving at the same time, at the same speed.

After positioning the animal on the related lane, the test proper is started by pushing the related **START** button/s (1 to 6 for mice, 1, 3 and 5 for rats), which start **the counters for each individual lane**.

The test proceeds until one of the conditions entered in **SETTING** is met (number of shocks, distance, time, see paragraph 5.4.1), or if <u>the animal rests on the activat-</u> ed electrified grid for more than 10 seconds.

The data are saved and will appear in the RESULT table.





Please note the - and + buttons at the bottom of the panel: these enable the user to increase or decrease the speed even while the test is running.


In constant speed mode, when the test in a specific lane is over, a fresh mouse can be positioned and a new test started.



5.5.2 Starting a Ramp

In the accelerating mode, when the **START** button is depressed, the Treadmill will pass from the IDLE state (motor and counters disabled), to the **initial speed**.





The test proper is started by the **START** button, which should be pressed as soon as the initial speed is reached.

In the RAMP mode, the individual lane buttons are not active, as the test starts simultaneously in all lanes (to guarantee that the whole ramp is exploited correctly).



The test proceeds until one of the conditions entered in the Setting is met (number of shocks, distance, time), or if the animal rests on the activated electrified grid for more than 10 seconds.

As soon as the FINAL SPEED is reached, the belt will keep rotating at constant speed and the timers will continue counting.

When the test is over, the related counters stop and the display shows the lane with a blue backdrop; the data are saved and will appear in the RESULT table.

The motor and the counters will stop by pressing the STOP button at any time during the experiment.

Press 🙂 to return to the Home page.



5.6 Results

Data are saved in the internal memory and they can be scrolled on the touch-screen. Reach the results section by depressing the "RESULTS" icon from the home page.



The home page is reached from any menu, but depressing the VV button



Monday, 20 July 2015 9:41	٢	ا 🖃 ي	3
DATE - TIME:		LANE:	
TREATMENT:	ID:		
PROTOCOL:	STAGE:	TRIAL:	
T-ELAPSED:	PERM:	SPEED:	
MODE:	DISTANCE:		
DATE - TIME:		LANE:	
TREATMENT:	ID:		
PROTOCOL:	STAGE:	TRIAL:	
T-ELAPSED:	PERM:	SPEED:	-
MODE:	DISTANCE:		

Scroll through the results by the UP \triangle and DOWN \checkmark arrows.

As in the example, all the Experiment-related information entered by the user appear near the blue headings

- DATE-TIME
- = LANE
- TREATMENT
- ID
- PROTOCOL
- STAGE
- TRIAL

The test results appear under the red headings:

- T-ELAPSED (Time)
- **PERM (Permanence)**
- SPEED
- MODE
- DISTANCE

Results can then be saved on the USB flash drive, and uploaded on the PC for further processing, see paragraph 5.8.1. When transferred to PC via USB drive, test results appear in extended version.

5.6.1 How to Export, Interpret and Erase the Experimental Data

During the experiments on the treadmill, at each contact of the animal with the corresponding shock grid, the device saves in the internal memory a set of data which characterize the recorded event.



The first set of columns reports the experiment tag, which enables a precise identification: **TREATMENT**, **PROTOCOL**, **STAGE**, **TRIAL**, **ID**, see paragraph 5.3. The **DATE** column supplies date and time (with a precision in the order of seconds) at which the event was recorded.

The following columns describe the events in details:

- LANE indicates the lane to which the record relates
- **T-ELAPSED** the column "TIME-ELAPSED" indicates after how long after the experiment start the shock even took place
- **PERM** the "SHOCK-PERMANENCE" column represents the time the animal spent on the shock grid before stepping back on the belt. If the animal stays for 10s or longer (cut-off time), the column will show 10s.

This is an "emergency" cut-off introduced to avoid damages to the animal.

- MODE indicates the mode active in the current experiment (Constant, Ramp, etc.)
- **SPEED** reports the speed at which the belt was running when the animal stepped on the grid (with our without shock)
- **DISTANCE** shows the distance walked by the animal from the beginning of the experiment to the moment when the animal stepped on the grid (with our without shock)
- **NUMSHOCK** this counter reports the total shocks the animal received since the beginning of the trial
- **INCLINATION** records the slope of the running belt during the test, if the information was entered when defining the experiment, <u>see paragraph 5.3</u>.

5.7 Important Note On Shock Frequency/Detection

In our model, the detection of the animal stepping on the grid at the bottom of the running surface, is carried out by the same circuit which administers the shock.

For this reason, the sampling rate is affected by the selected shock frequency (1, 2 or 3Hz).

To exemplify the above: if the freuency of 1Hz is set, the sampling rate will be 1 per second: this means that if the animal is quick enough to step over the grid twice within a second, it will be recorded as a single event. In any case, this should be a rare occurrence indeed and should not affect the correct event recording.



5.8 USB Storage



By depressing the USB icon, you will reach the USB Storage menu.



5.8.1 Save Data



By pressing this icon, it is possible to export the results on a USB-KEY.

Insert in the upper USB port the USB flash drive provided with standard package. Correct insertion is confirmed by the \bigcirc icon on the toolbar.

If you try to save the data without USB-KEY, the following dialog box will appear.

When the download procedure ends correctly, the following dialog box is displayed.

The file will be saved in .CSV format.



5.8.2 Load Protocol



When a USB flash drive, containing at least one protocol is connected to the device, the "LOAD PROTOCOL" button will open the following window, displaying all the loadable protocols found.

Use the arrows on the right of the screen to scroll the list up and down.



Having selected the desired protocol, if the protocol was properly loaded the confirmation box will show up, while an error box will appear if the file has been found corrupted.

After opening a protocol, the user will be able to test the animals following to the loaded protocol, combined with the information contained therein.



5.9 Utilities

The UTILITIES menu offers the following options:



5.9.1 Software Update

For software update you will receive the "TM_Update.hex" file by E-Mail or download it from our web site.

The file is to be saved on a USB key (FAT32 file system formatted), for example the standard Ugo Basile USB key provided with the Treadmill.

Proceed as follows:

IPDATE

- A. Copy the "TM_Update.hex" file on USB Key (file system FAT32);
- B. Switch ON the Treadmill
- **C.** Insert the USB KEY containing "TM_Update.hex" file into the UPPER USB port (see picture on the right)




- If the USB key is detected, the vice icon on the system bar will appear on the tool-bar
- If the system recognizes a new SW version, the icon will appear on the toolbar.



Only if the two icons are shown, it is possible to proceed with the update.

D. Select in sequence:





- E. Confirm by answering **OK** to the question "Update Treadmill application?
- F. Wait about 30 seconds for the device to reboot.





IMPORTANT:

do not remove the USB Key during the updating procedure!!!

- G. Wait about 30 seconds
- H. Remove the USB key from the upper USB port (<u>NB: do NOT remove the USB</u> <u>Hard disk from the lower USB port!</u>)
- I. Press the HOME button to check the new SW version.



5.9.2 Setting Date & Time

Select the "DATE-TIME" icon



Set date and time by using the the UP ▲ and DOWN ▼ arrows, then confirm by pressing OK.





5.9.3 Diagnostics



By selecting this function, the Treadmill will perform an self-diagnostics automatic test.



Then the self-diagnostics checks the motor performances, in accelerating mode:


The procedure can be stopped at any time, by pressing the SAVE REPORT TO USB button.

When the diagnostic sequence is completed, a report can be saved to USB



First of all, the machine tests the correct functioning of all the lane buttons.

Follow the instructions for all 6 buttons, then press **OK** .



5.9.4



This button resets the HD, which contains ALL THE EXPERIMENTAL DATA.

A message box will ask to confirm the action.

ATTENTION:

Reset HD

After confirming the procedure, it will be impossible to recover the data.



5.9.5 Cleaning



The cleaning option is very useful at the end of an experimental session to help the user cleaning the device. When selected it will start the motion of the carpet at a low speed (10 m/min) for 2 minutes, keeping the shock grid deactivated. A blade positioned below the running belt will remove the animal dropping. See also paragraph 6.3.


5.10 Signal Output

For some sophisticated application, and interface with other systems, it might be useful to obtain outputs for lane status, rotation and speed.

The 15-pin D-SUB connector (see picture and paragraph 3.8) provides different output trigger, in the form of TTL signals.

All pins, except pin 8, are TTL digital signals. Pin 8 is an analog output which provides the rotation speed.

Pins 14 and 15 are ground.



5.10.1 Pin-Out

Pin No.	DESCRIPTION	LEVEL LOW	LEVEL HIGH
Pin 1	LANE 1 SHOCK	NO SHOCK	SHOCK DELIVERED
Pin 2	LANE 2 SHOCK	NO SHOCK	SHOCK DELIVERED
Pin 3	LANE 3 SHOCK	NO SHOCK	SHOCK DELIVERED
Pin 4	LANE 4 SHOCK	NO SHOCK	SHOCK DELIVERED
Pin 5	LANE 5 SHOCK	NO SHOCK	SHOCK DELIVERED
Pin 6	LANE 6 SHOCK	NO SHOCK	SHOCK DELIVERED
Pin 7	EXPERIMENT STATUS	NOT RUNNING	RUNNING
Pin 8 (*)	ANALOG SPEED OUT	3m/min= 0.5V	100m/min = 10V
Pins 9 – 13	NOT USED		
Pins 14-15	GROUND		



IMPORTANT:

Please note the TTL output is not real-time. Maximum delay: 200ms.



6 MAINTENANCE

While any service of the instrument ought to be carried out by Ugo Basile personnel or by qualified personnel authorized by UGO BASILE organization, this manual section describes normal maintenance procedures which can be carried out at your facility.

<u>UNPLUG THE MAINS CORD BEFORE CARRYING OUT ANY</u> <u>MAINTENANCE JOB!</u>

6.1 Electrical

To inspect and/or replace the fuses, **disconnect the mains cable first!** Insert a miniature screwdriver in the slot indentation, see paragraph 3.5.2, and snap out the slide which houses the fuses. Use T5A fuses for operation at both 230 or 115 Volts. Snap in the fuse slide: the mechanical "*click*" ensures that it is locked.

6.2 Cleaning

The 47300 Treadmill was designed to make cleaning very easy. For general cleaning, **do not use organic solutions**, liable to impair the lane assembly, the acrylic components and touch-pad. Cotton wool and water can be used for cleaning purposes. For disinfection, use a non-alcoholic disinfectant, or H_2O_2 .

Each part needing specific cleaning can be conveniently disassembled. Please see the following paragraphs and follow the instruction carefully.

6.3 Running Belt

The running surface of the 47300 consists of a white alimentary-grade belt, easy-to-clean and soft enough to provide suitable grip, but not porous.

This specially selected material is easy to disinfect and requires minimal maintenance.

The cleaning feature, in the utility menu, see paragraph 5.9.5, provides a very efficient tool to be used at the end of each experimental session.



A blade positioned below the running belt will remove the animal droppings: enhance the action of the blade by adding a moderate quantity of disinfectant on the belt during the cleaning phase.

6.3.1 Aligning the Belt

Make sure that the belt is always centred to the chassis of the Treadmill. When the belt is misaligned with respect to the axis of the Treadmill, it is necessary to reset it by regulating the screws placed in the rear part of the Treadmill.

Here is a list of the steps to be followed to perform the aforementioned regulation:



1. Enter the START page and select the "RUN" button

2. Press the button "+" to increase the speed up to 20 m/min



- 3. Place yourself behind the Treadmill and establish whether the belt is misaligned towards left or right.
- 4. <u>If the belt drifts to the left side</u>, rotate the left regulation screw clockwise of a quarter of a turn and do the same for the screw on the right but counter-clockwise.
- 5. <u>If the belt drifts to the right side</u>, rotate the left regulation screw counterclockwise of a quarter of a turn and do the same for the screw on the right but clockwise. See following pictures:



6. Observe the behaviour of the belt for at least two minutes. If necessary, repeat the operations described at points 3, 4 and 5.

6.3.2 Regulating the Belt's Strain

It is recommended to regulate the strain of the belt in case it starts slipping on the rollers. This operation is very important to lengthen the lifetime of the components of the Treadmill.

During the regulation of the strain of the belt, keep the belt itself as slow as possible.



🔊 📑 🖥

1. Enter the START page and select the "RUN" button

2. Press the button "+" to increase the speed up to 20 m/min, to check the belt strain.

3. If the belt slips then <u>slow down the</u> <u>belt until it stops</u>, rotate both screws of a quarter of a turn in clockwise sense and, if necessary, repeat points 2 and 3.



SETUP

Monday, 20 July 2015

EXPERIMENT

9:41

4. Finally check if, after regulating the strain, the belt is correctly aligned. If not, please refer to paragraph 6.3.1.

6.3.3 Lane Assembly

The lane assembly can be easily removed for cleaning: just lift the whole assembly with due care, by holding it via the two recessed handles provided.



Figure 6 "Removing the Lane Assembly"





Do not use organic solutions, liable to impair the acrylic components. Cotton wool and water can be used for daily cleaning.

For disinfection, use a non-organic disinfectant, or H_2O_2 .

6.3.4 Drop Pan

The drop-pan positioned below the grid at the bottom of the running surface can be removed for cleaning.



<u>Please follow the steps in sequence</u>, to avoid any damage to the cleaning blade positioned between the belt and the drop-pan:

A. UNPLUG THE MAINS CORD

- **B.** First remove the lane assembly as shown in the picture on paragraph 6.3.1.
- C. Detach the cable connecting the grid/button assembly to the controller





D. Remove the whole grid/button assembly








Figure 7 "Removing the Drop-Pan"

After cleaning, reposition the drop-pan by repeating the same actions in reverse order: reassemble the grid, connect the cable, re-position the lane-assembly.

6.3.5 Grid Assembly

UNPLUG THE MAINS CORD BEFORE CLEANING THE GRID



To assure correct operation, shock delivery and detection, the grid needs to be clean.

When necessary, it can be easily removed, following the steps outlined below:

- A. First remove the lane assembly as shown in the picture on paragraph 6.3.1.
- **B.** Detach the cable connecting the grid/button assembly to the controller, then remove the whole grid/button assembly as explained in paragraph 6.3.4.
- C. Hold the grid assembly as shown in the picture, in order to keep the button array, and the electronic block away from water/liquids.
- D. Clean the grid thoroughly, with the help of an abrasive sponge, and a suitable disinfectant.
- E. Dry carefully
- F. Reassemble by repeating the above steps in reverse order.





Figure 8 "Cleaning the Grid-Assembly"

6.4 Lubrication

Our Treadmills <u>do not require</u> any lubrication: all its rotating components are lubricated for life.

6.5 Long Inactivity

The instrument does not require any particular maintenance after long inactivity, except cleaning.

6.6 Customer Support

For any further information you may desire concerning the use and/or maintenance of the device, please do not hesitate to contact our **service department** (or our local distributor) either directly of via our support page http://www.ugobasile.com/support.html :





UGO BASILE s.r.l. Via G. Di Vittorio 2 21036 GEMONIO – Varese, ITALY

Phone: +39 0332 744574



service@ugobasile.com logistics@ugobasile.com sales@ugobasile.com

Before sending any instrument to our factory for repair, please contact our logistics department to obtain a return authorization number (RMA) and shipping/packing instructions.

We may not be held responsible for damages during transport due to poor packing; whenever possible, please use the original packing.

7 47300 SPECIFICATIONS

General Commands 4"3 touch-screen (resistive) Read-out on the touch-screen **Power Requirement** Universal input 85-264 VAC, 50-60Hz, 40W max. Sound Level 80 dB ca., at maximum speed 10° to 40° C **Operating Temperature** Operation adjustable in the range 3 to 100m/min, in steps of 1m/min Speed Mode constant, ramp (accelerating), multi-step ramp (NEW!) Slope Manual tilting, positive (uphill) or negative (downhill) from -25° to +25°, in steps of 5° Shocker included from 0 to 2mA (in 0.1mA steps), 1, 2 or 3Hz Shock from the touch-screen Start/Stop Detection via incorporated electronic circuit Results Speed and distances (absolute and relative) Data Acquisition via dedicated X-PAD software, provided Data Portability by USB flash drive TTL Output Shock status (for each lane) and speed

Page 31 47300 Instruction Manual (Rev. 0)



Data	exported to Text, Excel or Pdf, or saved to cloud to Drop- Box, OneDrive, GoogleDrive
Configurations	exported to Text, Excel or Pdf, or saved to cloud to Drop- Box, OneDrive, GoogleDrive
Physical	
Total Weight	47302 & 47303 : 22Kg – 47300 : 27Kg
Shipping Weight	47302 & 47303 : 35Kg – 47300 : 40Kg approx.
Dimensions	56(w)x67(d)x35(h)cm
Packing Dimensions	47302 & 47303 : wooden crate, 77x65x63cm
	47300 : wooden crate, 82x71x57cm
Warranty	47300 is covered by a 24-month warranty

8 ORDERING INFORMATION

- 47300-001 Tapis-Roulant with touch-screen controller and shocker
- **47300-002** 3-Lane partition assembly for Rats, each lane 45x11x15(h)cm, with transparent cover
- E-AU 041 USB pen drive, including:

47300-302Instruction ManualX-PADSoftware Package

E-WP008 Mains Cord – Europe (or E-WP008-1 U.S.A. / E-WP008-2 U.K.)

Set of 2 fuses

47303	Mouse	Treadmill I	NG,	complete with following accessories
-------	-------	-------------	-----	-------------------------------------

- 47300-001 Tapis-Roulant with touch-screen controller and shocker
- **47300-003** 6-Lane partition assembly for Mice, each lane 45x5.5x7(h)cm, with transparent cover
- E-AU 041 USB pen drive, including:

47300-302Instruction ManualX-PADSoftware Package

E-WP008 Mains Cord – Europe (or E-WP008-1 U.S.A. / E-WP008-2 U.K.)

Set of 2 fuses



Working with both rats and mice? You should consider the Combination Package 47300, including both Mouse and Rat accessories:

47300	Combo Package: Mouse & Rat Treadmill NG , complete with follow- ing accessories		
47300-001	Tapis-Roulant with touch-screen controller and shocker		
47300-002	3-Lane partition assembly for Rats, each lane 45x11x15(h)cm, with transpar- ent cover		
47300-003	6-Lane partition assembly for Mice, each lane 45x5.5x7(h)cm, with transpar- ent cover		
E-AU 041	USB pen drive, including:		
	47300-302 X-PAD	Instruction Manual Software Package	
E-WP008	Mains Cord	I – Europe (or E-WP008-1 U.S.A. / E-WP008-2 U.K.)	

Set of 2 fuses

8.1 Optional

We offer a special model for tethered mice, with no lid, and higher walls:

- **47313 Mouse Treadmill NG**, complete with standard accessories and in alternative to 47300-003:
- **47300-013** 6-lane assembly (each lane 45x5.5, height 15cm, without lid, for tethered mice)

9 **BIBLIOGRAPHY** - Method Papers

- American Physiological Society: "Resource Book for the Design of Animal Exercise Protocols" February 2006
- O.J. Kemi et alia: "Intensity-Controlled Treadmill Running in Mice: Cardiacand Skeletal Muscle Hypertrophy" <u>J. Appl. Physiol</u>. 93: 1301-1309, 2002
- X.Q. Wang & G.W. Wang: "Effects of Treadmill Exercise Intensity on Spatial Working Memory and Long-Term Memory in Rats" Life Sciences 149: 96-103, 2016
- M. Shinozaki et alia: "Combined Treatment With Chondroitinase ABC and Treadmillrehabilitation for Chronic Severe Spinal Cord Injury in Adult Rats" <u>Neuroscience Res</u> 113: 37-47, 2016

Notes
\sim

Notes

r

Notes



CE CONFORMITY STATEMENT

Manufacturer	UGO BASILE srl
manufacturer	UCC DAGILL SIT

Address Via G. di Vittorio, 2 – 21036 Gemonio, VA, ITALY

Phone n. +39 0332 744574

Fax n. +39 0332 745488

We hereby declare that

Instrument. RODENT TREADMILL

Catalog number 47300 / 47302 / 47303

It is manufactured in compliance with the following European Union Directives and relevant harmonized standards

- 2014/35/UE relating to electrical equipment designed for use within certain voltage limits
- 2006/42/CE on machinery
- 2014/30/UE relating to electromagnetic compatibility
- 2011/65/UE 2015/863/UE on the restriction of the use of certain hazardous substances in electrical and electronic equipment

Account Manager	Adriano Basile
	Nome / Name
January 2017	AK
	Firmed Classical
Date	Firma 1 Signature
	1